



Azma Elixir Pajoo

Cat. No. AFTDX-1053

For Research Use Only

## Fast **DNrich** Blood Kit

### Introduction

The Fast **DNrich** Blood Kit provides all of the reagents necessary to extract DNA from whole blood. DNA purified with this kit is suitable for a variety of applications, including amplification and digestion with restriction endonucleases.

Fast <b>DNrich</b> Blood Kit components	Cat. No. AFBDX-1053-050
<b>TD Buffer</b>	<b>20 ml</b>
<b>Activator Reagent*</b>	<b>2 ml</b>
<b>VI Buffer</b>	<b>15 ml</b>
<b>SE Buffer</b>	<b>15 ml</b>
<b>Wash Buffer (conc.) *</b>	<b>15 ml</b>
<b>Elution Buffer</b>	<b>3 ml</b>
<b>column</b>	<b>50</b>
<b>Manual</b>	<b>1</b>

\* Refer to reminder, **Activator Reagent Preparation** and **Wash Buffer Preparation** before first use.

### Chemical Hazard

Always wear gloves and practice standard safety precautions while using the kit. Do NOT disinfect extraction waste in solutions containing bleach or any other form of acid. To clean any items contaminated with the reagent, simply soak in detergent and water to remove all traces of contamination before cleaning with bleach or acidic solutions.

### Reminder

Pre-set heater block at 65°C.

Prepare **Activator Reagent** immediately prior to use. Prepared Activator Reagent must be kept at 4°C.

Prepare **Wash Buffer** before first use.

### Activator Reagent Preparation

Add **2 ml** of molecular biology grade water to **Activator Reagent**, and vortex it well.

Note: Mark the check box on the bottle and write the date.

Note: For the best results, the prepared Activator Reagent should be used immediately. Prepared Activator Reagent can be stored for up to 3 months or 6 months at 4°C and -20°C, respectively.

### Wash Buffer Preparation

Add **35 ml** molecular biology grade **Absolute Ethanol** to **Wash Buffer** bottle before first use and mark the check box on it.



## PROTOCOL

### Step a: Sample Preparation

- a1 Transfer **100 µl** of homogenized **blood** sample into a 1.5 ml tube.

### Step b: Digestion & Lysis

- b1 Add **400 µl** of **TD Buffer** and **40 µl** of **Activator Reagent** to the sample tube and vortex vigorously.
- b2 Incubate at **65°C** for 30 to 45 minutes.  
**Note:** During incubation time, vortex the sample tube every 10 minutes.
- b3 Incubate at **85°C** for **10 minutes** and vortex every 5 minutes.
- b4 **(Optional)** Add **300 µl** of **VI Buffer** to sample tube and vortex vigorously and then keep at RT for **5 minutes**.
- b5 Centrifuge at **11500 g** for **10 minutes**.  
**Note:** Do not disturb the phases.
- b6 Carefully transfer about **300µl of supernatant** to a new 1.5 ml tube.
- b7 Add **300µl** of **SE Buffer** to the tube, invert for 5 times and keep at room temperature for **3 minutes** and then transfer all the sample to a **spin column**.
- b8 Centrifuge at **2000 g** for **2 minutes** and **discard** the flow through.

### Step c: Washing

- c1 Add **500 µl** of **Wash Buffer** to the column, centrifuge at **8000 g** for **2 minutes** and **discard** the flow through.
- c2 **Repeat step c1.**

### Step d: Column Drying

- d1 Centrifuge at **8000 g** for **1 minute**.
- d2 **Discard** the flow through and place the column into a new **1.5 ml** microcentrifuge tube.

### Step e: DNA Elution

- e1 Add **50 µl** of **Elution Buffer** to the center of column and let stay at **RT** for **3 minutes**.
- e2 Centrifuge at **10000 g** for **2 minutes**.

## Short PROTOCOL

### Step a: Sample Preparation

- a1 Transfer **100 µl** of homogenized **blood** sample into to a 1.5 ml tube.

### Step b: Digestion & Lysis

- b1 Add **400 µl** of **TD Buffer** and **40 µl** of **Activator Reagent** to the sample tube and vortex vigorously.
- b2 Incubate at **65°C** for 30 to 45 minutes.
- b3 Incubate at **85°C** for **10 minutes**.
- b4 Centrifuge at **11500 g** for **5 minutes**.
- b5 Carefully transfer about **300µl of supernatant** to a **spin column**.
- b6 Centrifuge at **2000 g** for **2 minutes** and **discard** the flow through.

### Step c: Washing

- c1 Add **500 µl** of **Wash Buffer** to the column, centrifuge at **8000 g** for **2 minutes** and **discard** the flow through.
- c2 **Repeat step c1.**

### Step d: Column Drying

- d1 Centrifuge at **8000 g** for **1 minute**.
- d2 **Discard** the flow through and place the column into a new **1.5 ml** microcentrifuge tube.

### Step e: DNA Elution

- e1 Add **50 µl** of **Elution Buffer** to the center of column and let stay at **RT** for **3 minutes**.
- e2 Centrifuge at **10000 g** for **2 minutes**.

